

Product Datasheet

Oligo (dT) Mag Beads (orb2663243)

Catalog Number	orb2663243
Category	Tools
Description	biorbyt's Oligo(dT) Magnetic Beads are designed for the rapid isolation of high-purity eukaryotic mRNA from total RNA or direct lysates via poly(A) tail hybridization. The purified mRNA is ideal for downstream applications including RT-PCR, cDNA library construction, and Northern blot analysis in both in vitro and in vivo research contexts.
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)