

Product Datasheet

Mouse IL-2 Antibody (orb2665700)

Catalog Number	orb2665700
Category	Antibodies
Description	This rat anti-mouse IL-2 monoclonal antibody (clone BU6D8X3, isotype IgG2a κ) is used for IL-2 neutralization and detection in both in vitro assays and in vivo mouse models of immunology and T-cell biology research. The recommended isotype control is Rat IgG2a κ .
Target	IL Receptor
Clonality	Monoclonal
Conjugation	Unconjugated
Concentration	8.61 mg/mL
Purity	95.00% (May vary between batches)
MW	150 kDa
Antibody Type	Inhibitor Antibody
Clone Number	BU6D8X3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)