

Product Datasheet

CEACAM7 Antibody (orb2665709)

Catalog Number	orb2665709
Category	Antibodies
Description	BR5X9D is a mouse IgG1 kappa chimeric monoclonal antibody targeting human CEACAM7, supplied with a mouse IgG1 kappa isotype control. It is a research tool used in vitro for the detection and functional study of CEACAM7 in colorectal cancer and other epithelial cell research.
Target	Others
Clonality	Monoclonal
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Clone Number	BR5X9D
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)