

Product Datasheet

Vadastuximab Antibody (orb2665784)

Catalog Number	orb2665784
Description	Vadastuximab (Anti-Human CD33 Recombinant Antibody) is a humanized anti-CD33 monoclonal antibody that can be used for the synthesis of the ADC compound Vadastuximab talirine and for the study of acute myeloid leukemia (AML).
Conjugation	Unconjugated
Target	Vadastuximab
Concentration	6.32 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	98.00%
MW	145.5 kDa
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)