

Product Datasheet

PE-conjugated Anti-CD40 Ligand antibody(DM164); Rabbit mAb (orb2669184)

Catalog Number	orb2669184
Description	PE-conjugated Anti-CD40 Ligand antibody(DM164); Rabbit mAb
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Target	CD40 Ligand
Storage	Store at 2°C-8°C for 6 months
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Clone Number	DM164
Uniprot ID	P29965
Dilution Range	Flow Cyt 1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)