

Product Datasheet

Canine IL31RA Biosimilar Antibody (orb2669629)

Catalog Number	orb2669629
Category	Antibodies
Description	Anti-Canine IL31RA(nolavetbart biosimilar) IgG2) mAb
Target	Canine IL31RA
Clonality	Monoclonal
Isotype	IgG2
Conjugation	Unconjugated
Reactivity	Canine
Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
UniProt ID	XP_013963900.1
Tested applications	FC
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558