

## Product Datasheet

### LILRB4 Biosimilar Antibody (orb2669889)

<b>Catalog Number</b>	orb2669889
<b>Category</b>	Antibodies
<b>Description</b>	Biotinylated Anti-LILRB4(rolistobart biosimilar) mAb
<b>Target</b>	LILRB4
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG4
<b>Conjugation</b>	Biotin
<b>Reactivity</b>	Human
<b>Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>UniProt ID</b>	<b>Q8NHJ6</b>
<b>Tested applications</b>	FC
<b>Antibody Type</b>	Biosimilar Antibody
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)