

Product Datasheet

TSLP Biosimilar Antibody (orb2670072)

| | |
|----------------------------|--|
| Catalog Number | orb2670072 |
| Category | Antibodies |
| Description | Biotinylated Anti-TSLP(ocankitug biosimilar) mAb |
| Target | TSLP |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Conjugation | Biotin |
| Reactivity | Human |
| Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| UniProt ID | Q969D9 |
| Tested applications | FC |
| Antibody Type | Biosimilar Antibody |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)