

Product Datasheet

ITGA4 Biosimilar Antibody (orb2670818)

Catalog Number	orb2670818
Category	Antibodies
Description	PE-conjugated Anti-ITGA4 (natalizumab biosimilar) mAb
Target	ITGA4
Clonality	Monoclonal
Species/Host	Human Humanized
Isotype	IgG4
Conjugation	Unconjugated
Reactivity	Human
UniProt ID	P13612
Tested applications	FC
Dilution range	Flow Cyt 1:100
Antibody Type	Biosimilar Antibody
Storage	Store at 2°C-8°C for 6 months
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)