

Product Datasheet

BAFF-R Biosimilar Antibody (orb2670841)

Catalog Number	orb2670841
Category	Antibodies
Description	PE-conjugated Anti-BAFFR (ianalumab biosimilar) mAb
Target	BAFF-R
Clonality	Monoclonal
Species/Host	Human
Isotype	IgG1
Conjugation	PE
Reactivity	Human
Reconstitution	Liquid □ PBS with 0.18% BSA □ 0.1% Proclin300
UniProt ID	Q96RJ3
Tested applications	FC
Dilution range	1 μl/test
Antibody Type	Biosimilar Antibody
Storage	Store at 2°C-8°C for 6 months
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)