

Product Datasheet

GABARAPL2 Antibody (orb2673417)

Catalog Number	orb2673417
Category	Antibodies
Description	GABARAPL2 Antibody
Target	GABARAPL2
Clonality	Polyclonal
Conjugation	Unconjugated
Reactivity	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to a.a. 1-117 of human GABARAPL2 (NP_009216.1).
Tested applications	IF, IHC, WB
Dilution range	WB 1:500 - 1:2000, IHC 1:50 - 1:200 IF 1:20 - 1:50
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)