

Product Datasheet

RCE1 Antibody (orb2677571)

Catalog Number	orb2677571
Category	Antibodies
Description	RCE1 Antibody
Target	RCE1
Clonality	Polyclonal
Conjugation	Unconjugated
Reactivity	Human
Immunogen	Recombinant protein of human RCE1.
Tested applications	IF, IHC, WB
Dilution range	WB 1:500 - 1:2000, IHC 1:50 - 1:200 IF 1:50 - 1:200
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)