

Product Datasheet

UBR5 Antibody (orb2681979)

Catalog Number	orb2681979
Category	Antibodies
Description	UBR5 Antibody
Target	UBR5
Clonality	Polyclonal
Conjugation	Unconjugated
Reactivity	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to a.a. 550-760 of human UBR5 (NP_056986.2).
Tested applications	WB
Dilution range	WB 1:500 - 1:2000
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)