

Product Datasheet

TIP 39 (39 mer) (orb2692903)

Catalog Number	orb2692903
Category	Proteins
Description	TIP 39, Tuberoinfundibular Neuropeptide is a neuropeptide and parathyroid hormone 2 receptor (PTH2R) agonist. TIP 39 is highly conserved among species. TIP39 from all species activates adenylyl cyclase and elevates intracellular calcium levels through parathyroid hormone 2 receptor (PTH2R).
Target	Adenylate Cyclase
Purity	≥95%
Protein Sequence	SLALADDAAFRERARLLAALERRHWLNSYMHKLLVLDAP
MW	4504.28
Solubility (25°C)	Water: 50 mg/mL
CAS Number	277302-47-3
Formula	C ₂₀₂ H ₃₂₅ N ₆₁ O ₅₄ S
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)