

Product Datasheet

(Pro3) Gastric Inhibitory Peptide (GIP), human (orb2693935)

Catalog Number	orb2693935
Category	Proteins
Description	(Pro3) GIP, human ((Pro3) Gastric Inhibitory Peptide, human) is an efficacious, stable and specific human GIP receptor (hGIPR) full agonist. (Pro3) GIP, human has high binding affinity for human GIPR with Ki/ Kd values of 0.90 nM. (Pro3) GIP, human can be used for the research of obesity-related diabetes.
Target	Insulin Receptor
Purity	≥95%
Protein Sequence	YAPGTFISDYSIAMDKIHQQDFVNWLLAQKGKKNDWKHNITQ
MW	4951.53
CAS Number	299898-52-5
Formula	$C_{226}H_{338}N_{60}O_{64}S$
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)