

## Product Datasheet

### PAR-1-selective peptide (orb2693980)

<b>Catalog Number</b>	orb2693980
<b>Category</b>	Proteins
<b>Description</b>	Protease-Activated Receptor-1, PAR-1 Agonist is a selective proteinase-activated receptor1 (PAR-1) agonist peptide. Protease-Activated Receptor-1, PAR-1 Agonist corresponds to PAR1 tethered ligand and which can selectively mimic the actions of thrombin via this receptor.
<b>Target</b>	Protease Activated Receptor (PAR)
<b>Purity</b>	≥95%
<b>Protein Sequence</b>	TFLLRN
<b>MW</b>	762.91
<b>CAS Number</b>	141136-85-8
<b>Formula</b>	$C_{35}H_{58}N_{10}O_9$
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)