

Product Datasheet

Human IL12A & IL27B Protein, His & Fc Tag (orb2698055)

Description	Human IL12A & IL27B Protein, His & Fc Tag
Reactivity	Human
Endotoxins	> (59.8+33.3)%, determined by SDS-PAGE
Target	IL12A & IL27B
Note	For research use only
Purity	> (59.8+33.3)%, determined by SDS-PAGE
Source	<p>This protein contains the p35 subunit of human IL12, termed as IL12A (NP_0008732) (Met1-Ser219) was fused with the C-terminal His-tagged Fc region of human IgG1 at the C-terminus, constructed the plasmid 1; This protein contains the human IL27B, (NP_0057462) (Met1-Lys229) was fused with the C-terminal flag-tagged Fc region of human IgG1 at the C-terminus, constructed the plasmid 2 . The two plasmids were co-expressed and the heterodimer was purified and expressed from HEK293 Cells.</p>
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)