

Product Datasheet

Human IL12A & IL12B Protein, Fc Tag (orb2698078)

Description	Human IL12A & IL12B Protein, Fc Tag
Reactivity	Human
Endotoxins	> 95%, determined by SDS-PAGE
Target	IL12A & IL12B
Note	For research use only
Purity	> 95%, determined by SDS-PAGE
Source	This protein contains the p35 subunit of human IL-12, termed as IL12A (NP_000873.2) (Met 1-Ser 219) was fused with a polyhistidine tag at the C-terminus, constructed the plasmid 1 ; This protein contains the p40 subunit of human IL-12, termed as IL12B (NP_002178.2) (Met 1-Ser 328) was fused with a polyhistidine tag at the C-terminus, constructed the plasmid 2. The two plasmids were co-expressed and the heterodimer was purified.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)