

Product Datasheet

Human SP-D/Surfactant protein D ELISA Kit (DIY Antibody Pairs) (orb2702476)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number | orb2702476 |
| Category | Assays and Kits |
| Description | Human SP-D/Surfactant protein D ELISA Kit (DIY Antibody Pairs) |
| Target | Pulmonary surfactant-associated protein D |
| Reactivity | Human |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: A21-F375 |
| UniProt ID | P35247 |
| MW | 37728 Da |
| Tested applications | ELISA |
| Specificity | Natural and recombinant human SP-D |
| Cross Reactivity | There is no detectable cross-reactivity with other relevant proteins. |
| Range | 625 pg/ml - 40,000 pg/ml |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)