

Product Datasheet

Human CCL24/Eotaxin-2 ELISA Kit (DIY Antibody Pairs) (orb2702479)

Catalog Number	orb2702479
Category	Assays and Kits
Description	Human CCL24/Eotaxin-2 ELISA Kit (DIY Antibody Pairs)
Target	C-C motif chemokine 24
Reactivity	Human
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V27-C119
UniProt ID	O00175
MW	13134 Da
Tested applications	ELISA
Specificity	Natural and recombinant human CCL24
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	31.2 pg/ml - 2,000 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)