

## Product Datasheet

# Human IL-1 RL1/ST2 ELISA Kit (DIY Antibody Pairs) (orb2702547)

<b>Catalog Number</b>	orb2702547
<b>Category</b>	Assays and Kits
<b>Description</b>	Human IL-1 RL1/ST2 ELISA Kit (DIY Antibody Pairs)
<b>Target</b>	Interleukin-1 receptor-like 1
<b>Reactivity</b>	Human
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: K19-S328
<b>UniProt ID</b>	<b>Q01638</b>
<b>MW</b>	63358 Da
<b>Tested applications</b>	ELISA
<b>Specificity</b>	Natural and recombinant human IL1RL1 as well as free IL1RL1 and IL-33 complexed IL1RL1.
<b>Cross Reactivity</b>	There is no detectable cross-reactivity with other relevant proteins.
<b>Range</b>	156 pg/ml - 10,000 pg/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)