

Product Datasheet

Human RBP4/Retinol Binding Protein 4 ELISA Kit (DIY Antibody Pairs) (orb2702552)

Catalog Number	orb2702552
Category	Assays and Kits
Description	Human RBP4/Retinol Binding Protein 4 ELISA Kit (DIY Antibody Pairs)
Target	Retinol-binding protein 4
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: E19-L201
UniProt ID	P02753
Tested applications	ELISA
Specificity	Natural and recombinant human RBP4
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	312 pg/ml - 20,000 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum, plasma(heparin, EDTA) and urine.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)