

Product Datasheet

DEDN6526A Biosimilar (ETBR Reference Antibody) (orb2703459)

Catalog Number	orb2703459
Category	Antibodies
Description	DEDN6526A (RG-7636) is a humanized antibody expressed in CHO cells, targeting ETBR. DEDN6526A (RG-7636) has a hulgG1 type heavy chain and a huk type light chain, with a predicted molecular weight (MW) of 145 kDa. The isotype control for DEDN6526A (RG-7636) can be referenced as Human IgG1 kappa, Isotype Control (HY-P99001).
Target	DEDN6526A Biosimilar(Anti-ETBR Reference Antibody)
Clonality	Monoclonal
Conjugation	Unconjugated
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)