

Product Datasheet

Ecromeximab Biosimilar (GD3 Reference Antibody) (orb2703715)

Catalog Number	orb2703715
Category	Antibodies
Description	Ecromeximab (KM 871; KW 2871) is a mouse-derived IgG-κ chimeric antibody targeting ganglioside GD3. Ecromeximab is expressed by YB2/0 rat hybridoma cells.
Target	Ecromeximab Biosimilar □ Anti-GD3 Reference Antibody □
Clonality	Monoclonal
Conjugation	Unconjugated
Antibody Type	Biosimilar Antibody
CAS Number	292819-64-8
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)