

## Product Datasheet

# Ficlatuzumab Biosimilar (Anti-HGF / SF Reference Antibody) (orb2703783)

|                       |                                                                                                                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b> | orb2703783                                                                                                                                                                                                                                                   |
| <b>Category</b>       | Antibodies                                                                                                                                                                                                                                                   |
| <b>Description</b>    | Ficlatuzumab is a monoclonal antibody (McAb) targeting human hepatocyte growth factor (HGF). Ficlatuzumab blocks the activation of the HGF/c-Met signaling pathway, and inhibits c-Met receptor-mediated cancer cell proliferation, migration, and invasion. |
| <b>Target</b>         | Ficlatuzumab BiosimilarAnti-HGF / SF Reference Antibody                                                                                                                                                                                                      |
| <b>Clonality</b>      | Monoclonal                                                                                                                                                                                                                                                   |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                                                                                                                                 |
| <b>Antibody Type</b>  | Biosimilar Antibody                                                                                                                                                                                                                                          |
| <b>CAS Number</b>     | 1174900-84-5                                                                                                                                                                                                                                                 |
| <b>Storage</b>        | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                                                                                                      |
| <b>Note</b>           | For research use only                                                                                                                                                                                                                                        |

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)