

Product Datasheet

SAIT301 Biosimilar (HGFR / c-Met Reference Antibody) (orb2703789)

Catalog Number	orb2703789
Category	Antibodies
Description	SAIT301 Biosimilar (HGFR / c-Met Reference Antibody)
Target	SAIT301 Biosimilar Anti-HGFR / c-Met Reference Antibody
Clonality	Monoclonal
Conjugation	Unconjugated
Purity	>98%
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558