

Product Datasheet

Anifrolumab Biosimilar (IFNAR1 Reference Antibody) (orb2703811)

| | |
|------------------------|--|
| Catalog Number | orb2703811 |
| Category | Antibodies |
| Description | Anifrolumab is a type I interferon (IFN) receptor antagonist, a human monoclonal antibody. Anifrolumab blocks the activity of type I interferon. Anifrolumab can be used in systemic lupus erythematosus (SLE) research. |
| Target | Anifrolumab Biosimilar(Anti-IFNAR1 Reference Antibody) |
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purity | >98% |
| Antibody Type | Biosimilar Antibody |
| CAS Number | 1326232-46-5 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)