

Product Datasheet

Astegolimab Biosimilar (IL-1RL1 / ST2 / IL-33R Reference Antibody (orb2703855)

Catalog Number	orb2703855
Category	Antibodies
Description	Astegolimab (MSTT 1041A; RG 6149) is a human IgG2 monoclonal antibody that blocks IL-33 signaling by targeting ST2, the IL-33 receptor. Astegolimab has the potential for chronic obstructive pulmonary disease (COPD) research.
Target	Astegolimab Biosimilar (Anti-IL-1RL1 / ST2 / IL-33R Reference Antibody)
Clonality	Monoclonal
Conjugation	Unconjugated
Purity	>98%
Antibody Type	Biosimilar Antibody
CAS Number	2173054-79-8
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)