

Product Datasheet

Mibavademab Biosimilar (Anti-LEPR / CD295 Reference Antibody) (orb2703963)

Catalog Number	orb2703963
Category	Antibodies
Description	Mibavademab Biosimilar (Anti-LEPR / CD295 Reference Antibody)
Target	Mibavademab Biosimilar(Anti-LEPR / CD295 Reference Antibody)
Clonality	Monoclonal
Conjugation	Unconjugated
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)