

## Product Datasheet

# RG7841 Biosimilar (Ly6E Reference Antibody) (orb2703972)

<b>Catalog Number</b>	orb2703972
<b>Category</b>	Antibodies
<b>Description</b>	RG7841 Biosimilar (Ly6E Reference Antibody)
<b>Target</b>	RG7841 Biosimilar Anti-Ly6E Reference Antibody
<b>Clonality</b>	Monoclonal
<b>Conjugation</b>	Unconjugated
<b>Purity</b>	>98%
<b>Antibody Type</b>	Biosimilar Antibody
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558