

## Product Datasheet

# Zelminemab Biosimilar (PAC1 Reference Antibody) (orb2704032)

<b>Catalog Number</b>	orb2704032
<b>Category</b>	Antibodies
<b>Description</b>	Zelminemab (AMG-301) is a humanized monoclonal antibody that inhibits PACAP type I (PAC1) receptor.
<b>Target</b>	Zelminemab Biosimilar(Anti-PAC1 Reference Antibody)
<b>Clonality</b>	Monoclonal
<b>Conjugation</b>	Unconjugated
<b>Purity</b>	>98%
<b>Antibody Type</b>	Biosimilar Antibody
<b>CAS Number</b>	2225850-33-7
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)