

Product Datasheet

Phospho-MCM2 (Ser27) Rabbit mAb (orb2707003)

| | |
|-----------------------------|---|
| Catalog Number | orb2707003 |
| Category | Antibodies |
| Description | Phospho-MCM2 (Ser27) Rabbit mAb |
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Reactivity | Human, Mouse, Rat |
| Form/Appearance | Liquid |
| Buffer/Preservatives | 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% rAlbumin, 0.02% sodium azide and 50% glycerol. |
| UniProt ID | P49736 |
| Tested applications | ICC, IF, IHC-P, IP, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)