

Product Datasheet

Rabbit Anti-Sheep IgG H&L(APC) (orb2714427)

Catalog Number	orb2714427
Category	Antibodies
Description	Rabbit Anti-Sheep IgG H&L(APC)
Species/Host	Rabbit
Conjugation	APC
Reactivity	Sheep
Form/Appearance	Liquid
Buffer/Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Tested applications	ELISA, FC, ICC, IF, IHC-Fr, IHC-P
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)