

Product Datasheet

Rabbit Anti-Sheep IgG H&L (Biotin) (orb2714450)

Catalog Number	orb2714450
Category	Antibodies
Description	Rabbit Anti-Sheep IgG H&L (Biotin)
Species/Host	Rabbit
Conjugation	Biotin
Reactivity	Sheep
Form/Appearance	Liquid
Buffer/Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Tested applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)