

## Product Datasheet

# Rabbit Goat IgG (H+L) Secondary Antibody (Cy3) (orb27706)

<b>Catalog Number</b>	orb27706
<b>Category</b>	Antibodies
<b>Description</b>	Rabbit Goat IgG (H+L) Secondary Antibody (Cy3)
<b>Clonality</b>	Polyclonal
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Cy3
<b>Reactivity</b>	Goat
<b>Concentration</b>	1mg/ml
<b>Tested applications</b>	FC, IF
<b>Specificity</b>	This Cy3 conjugated antibody is specific for goat IgG and shows no cross-reactivity with human/rat/mouse/rabbit IgG.
<b>Antibody Type</b>	Secondary Antibody
<b>Storage</b>	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)