

Product Datasheet

IL-23R Protein, Human, Recombinant (His & Avi), Biotinylated (orb2802696)

Catalog Number	orb2802696
Description	IL-23R Protein, Human, Recombinant (His & Avi), Biotinylated is expressed in HEK293 Cells. The accession number is Q5VWK5.
Storage	-20°C
Tag	C-hFc-Avi
Note	For research use only
Application notes	Please refer to the lot-specific COA.
MW	66.41 kDa (predicted); 95.6 kDa (reducing contition)
Uniprot ID	Q5VWK5
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)