

Product Datasheet

Pasritamig Antibody (orb2803424)

Catalog Number	orb2803424
Category	Antibodies
Description	Pasritamig is a bispecific T-cell engager antibody designed to target kallikrein 2 (KLK2) and CD3E. Its IgG1-scFv format facilitates the redirection of T-cells against KLK2-expressing cells, supporting research in prostate cancer immunotherapy for in vitro cytotoxicity assays and in vivo efficacy studies.
Target	Others
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)