

Product Datasheet

Mipletamig Antibody (orb2803508)

Catalog Number	orb2803508
Category	Antibodies
Description	Mipletamig is a bispecific chimeric antibody targeting CD3E, constructed from human and humanized components with a [(scFv - κ -heavy) G1-h-CH2-CH3 (scFv-heavy - κ)] ₂ format. It is a research tool used in vitro and in vivo to study T-cell engagement and activation in immunology and oncology research.
Target	Integrin
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)