

Product Datasheet

Avitotamig Antibody (orb2803519)

Catalog Number	orb2803519
Category	Antibodies
Description	Avitotamig is a bispecific T-cell engager antibody targeting CD3 and CD33. Its [H-γ1 _ L-κ-scFv heavy -κ]-dimer format facilitates T-cell-mediated cytotoxicity against CD33-expressing cells, supporting research in acute myeloid leukemia for in vitro and in vivo experimental models.
Target	Integrin
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)