

Product Datasheet

Spevatamig Antibody (orb2803540)

Catalog Number	orb2803540
Category	Antibodies
Description	Spevatamig is a bispecific antibody with a human IgG1/kappa format that simultaneously targets CD47 and CLDN18. It is utilized in research to investigate potential therapeutic strategies in oncology, particularly in preclinical in vitro and in vivo models for solid tumors.
Target	Immunology/Inflammation related
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)