

Product Datasheet

Androgen Receptor Antibody (orb2803569)

Catalog Number	orb2803569
Category	Antibodies
Description	Androgen Receptor Antibody
Target	Androgen Receptor
Clonality	Recombinant
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human, Mouse
Concentration	1000 µg/ml
Buffer/Preservatives	Phosphate Buffered Saline (PBS) with 0.09% Sodium Azide, rAlbumin-Free
Immunogen	between 1 and 50
UniProt ID	P10275
Tested applications	FC, ICC, IHC, IP, WB
Application notes	Format: Whole IgG
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)