

Product Datasheet

Human IL-1 Alpha/IL-1F1/IL1A Quick ELISA Kit (orb2806267)

Catalog Number	orb2806267
Category	Assays and Kits
Description	Human IL-1 Alpha/IL-1F1/IL1A Quick ELISA Kit
Reactivity	Human
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S113-A271
UniProt ID	P01583
MW	31 kDa
Tested applications	ELISA
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	4.7 pg/ml - 300 pg/ml
Sensitivity	0.5 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Storage	Store at 4°C for 6 months. (Ships with gel ice, can store for up to 3 days in room temperature. Refrigerate upon receipt.)
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)