

## Product Datasheet

### MST4 Antibody (orb2808361)

|                      |                                                                                              |
|----------------------|----------------------------------------------------------------------------------------------|
| Catalog Number       | orb2808361                                                                                   |
| Category             | Antibodies                                                                                   |
| Description          | MST4 Antibody                                                                                |
| Target               | MST4                                                                                         |
| Clonality            | Monoclonal                                                                                   |
| Species/Host         | Rabbit                                                                                       |
| Isotype              | Rabbit IgG                                                                                   |
| Conjugation          | Unconjugated                                                                                 |
| Reactivity           | Human, Mouse, Rat                                                                            |
| Concentration        | Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration. |
| Buffer/Preservatives | Purified antibody in PBS with 0.05% sodium azide, 0.05% rAlbumin and 50% glycerol.           |
| Immunogen            | A synthesized peptide derived from human MST4                                                |
| UniProt ID           | <b>Q9P289</b>                                                                                |
| MW                   | Calculated MW: 47 kDa ; Observed MW: 50 kDa                                                  |
| Tested applications  | IHC, WB                                                                                      |
| Dilution range       | WB: 1/1000-1/2000; IHC: 1/100-1/200                                                          |
| Antibody Type        | Primary Antibody                                                                             |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)