

Product Datasheet

Human 8-Epi Prostaglandin F2 Alpha (8-epi-PGF2a) ELISA Kit (orb2809663)

Catalog Number	orb2809663
Category	Assays and Kits
Description	Human 8-Epi Prostaglandin F2 Alpha (8-epi-PGF2a) ELISA Kit
Reactivity	Human
Application notes	Assay Type: Competitive
Assay Type	Competitive
Range	31.2-2000pg/mL
Sensitivity	15.1pg/mL
Sample Types	Serum;plasma and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)