

Product Datasheet

Human Heat shock-related 70 kDa Protein 2 (HSPA2) ELISA Kit (orb2809827)

Catalog Number	orb2809827
Category	Assays and Kits
Description	Human Heat shock-related 70 kDa Protein 2 (HSPA2) ELISA Kit
Reactivity	Human
UniProt ID	P54652
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	1.56-100ng/mL
Sensitivity	0.75ng/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)