

Product Datasheet

Mouse Dopamine (DA) ELISA Kit (orb2809930)

Catalog Number	orb2809930
Category	Assays and Kits
Description	Mouse Dopamine (DA) ELISA Kit
Reactivity	Mouse
Application notes	Assay Type: Competitive
Assay Type	Competitive
Range	15.6-1000pg/mL
Sensitivity	7.4pg/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)