

Product Datasheet

Mouse Phosphorylation Adenosine Monophosphate Activated Protein Kinase (p-AMPK) ELISA Kit (orb2810119)

| | |
|--------------------------|--|
| Catalog Number | orb2810119 |
| Category | Assays and Kits |
| Description | Mouse Phosphorylation Adenosine Monophosphate Activated Protein Kinase (p-AMPK) ELISA Kit |
| Reactivity | Mouse |
| Application notes | Assay Type: Sandwich |
| Assay Type | Sandwich |
| Range | 1.56-100ng/mL |
| Sensitivity | 0.72ng/mL |
| Sample Types | Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)