

Product Datasheet

Mouse Adenosine Monophosphate Activated Protein Kinase (AMPK) ELISA Kit (orb2810120)

Catalog Number	orb2810120
Category	Assays and Kits
Description	Mouse Adenosine Monophosphate Activated Protein Kinase (AMPK) ELISA Kit
Reactivity	Mouse
UniProt ID	Q9R078
Application notes	Assay Type: Sandwich
Assay Type	Sandwich
Range	1.56-100ng/mL
Sensitivity	0.74ng/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)