

Product Datasheet

Human Epinephrine (EPI) ELISA Kit (orb2810257)

Catalog Number	orb2810257
Category	Assays and Kits
Description	Human Epinephrine (EPI) ELISA Kit
Reactivity	Human
Application notes	Assay Type: Competitive
Assay Type	Competitive
Range	31.2-2000pg/mL
Sensitivity	15.4pg/mL
Sample Types	Serum;plasma;tissue homogenates;cell lysates;cell culture supernates and other biological fluids
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)