

Product Datasheet

Human Peroxisome proliferator-activated receptor gamma coactivator 1-alpha (PGC-1 α) ELISA Kit (orb2810569)

Catalog Number	orb2810569
Description	Human Peroxisome proliferator-activated receptor gamma coactivator 1-alpha (PGC-1 α) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	Assay Type: Sandwich
Uniprot ID	Q9UBK2
Sensitivity	0.08ng/mL
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)